

TREND

TREND DCT-QF37

QUICKFORM LSO

FEATURES:

- The **Trend Pro Quick Form** is a one-size, extremely easy to adjust LSO.
- **Patient Friendly.** This brace is the ideal product for your patients who have trouble adjusting products themselves.
- We have designed this product to go from **25"-60" with minimal effort**. This brace is very light weight, breathable, and offers excellent support and compression.



MOISTURE WICKING

The brace features moisture wicking material and breathable zones for a more comfortable patient experience.



EFFORTLESS SIZING

Easy adjust sizing guides allow for quick sizing and patient to apply correct pressure every time.



PDAC APPROVED
L0637/ L0650

NEW PULLEY SYSTEM

Revolutionary SYSTEM COMPRESSION

WHAT IS IT?

The Trend LSO offers a revolutionary design and compression system to easily mold and form to patients back. It encompasses a dual or single pull draw system that is great for patients that lack strength in both arms. The design offers patients a comfortable and easy to use LSO with superior quality, comfort and performance.



Quick Sizing Guide Sizing System

HOW DOES IT WORK?

By removing the waist wings from their velcro attachments, it reveals a simple sizing system that allows easy readjustment. Place on the estimated size, then replicate the process on the other side of the brace.



Contact Us For More Information.

(800) 441-2591

TREND-MED.COM